

CLEARER VISUALIZATION. SIMPLIFIED DISSECTION.¹



ProdiGI™ traction wire.
Bringing confidence to Endoscopic
Submucosal Dissection (ESD)

Medtronic

DESIGNED WITH YOU IN MIND



The ProdiGI™ traction wire has been purpose-built for ease-of-use — to deliver uninterrupted visualization of the dissection plane — bringing clarity and confidence to your ESD procedures.¹ It allows you to hold tissue during dissection procedures and expose the dissection plane.

ESD Simplified

In an ex-vivo study comparing against conventional ESD, ProdiGI™ traction wire:

- Significantly improved endoscopists' visualization of the submucosal space¹
- Exceeded physician expectations for ease of placement¹
- Significantly improved physician satisfaction with respect to time and traction effect¹
- Successfully used in esophageal, gastric and colorectal ESD procedures¹

Indicated Use

The Medtronic ProdiGI™ Traction Wire is indicated to grasp tissue within the gastrointestinal tract during an Endoscopic Submucosal Dissection procedure.

Complications

Potential complications include:

- Irritation or injury to the mucosa
- Device fragments in patient
- Tissue damage
- Allergic reaction
- Fever
- Mucosal laceration
- Minor acute bleeding
- Major bleeding
- Transfusion secondary to major bleeding
- Infection
- Perforation and secondary complications related to perforation
- Surgery or secondary intervention to correct perforation
- Aspiration
- Peritonitis
- Tissue seeding
- Death
- Additional complications associated with an endoscopic submucosa dissection procedure

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.06/2020-US-DG-2000029

5920 Longbow Drive
Boulder, CO 80534

800.255.8522
303.581.6898

medtronic.com/gi

Specifications

Model No	Device working lengths	Traction wire length	Endoscope working channel diameter
ERD-TW20	230 cm	20 mm	2.8 mm or larger
ERD-TW35	230 cm	35 mm	2.8 mm or larger
ERD-TWSG	230 cm	NA	2.8 mm or larger

Item No

Order code	Product description	Selling unit	UDI #	Box dimensions
ERD-TW20	ProdiGI™ traction wire 20 mm	Box of 1 system	10884521806863	12.5" x 11.3" x 1.3"
ERD-TW35	ProdiGI™ traction wire 35 mm	Box of 1 system	10884521806870	12.5" x 11.3" x 1.3"
ERD-TWSG	ProdiGI™ secondary grasper	Box of 2 graspers	20884521808383	12.5" x 11.3" x 1.3"

Ordering Information

Vendor:
Covidien Sales, LLC
3555 Koger Blvd Suite 200
Duluth, GA 30096

Phone: (770) 662-0870
Fax: (770) 662-0510

Email to send electronic PO:
gi-usacustomer@covidien.com

TAX ID: 02-0502162

References

1. RE00221182 (A) ProdiGI™ Traction Wire Design Validation Report.

Medtronic